
CYCLOHEXENE



Cyclohexene concentration (ppm)



1. PERFORMANCE

7) Colour change

1) Measuring range : 20-300 ppm Number of pump strokes : 1(100mL)

2) Sampling time : 1.5 minutes / 1 pump stroke

3) Detectable limit
4) Shelf life
5) Operating temperature
2 years
15 ~25℃
15 ~25℃

6) Reading : The tube scales are calibrated based on Ethyl acetate at 1 pump stroke and

Cyclohexene concentration is determined by using a conversion chart

at 1 pump stroke ∴ Yellow → Brown

2. RELATIVE STANDARD DEVIATION

RSD-low: - RSD-mid.: - RSD-high: -

3. CHEMICAL REACTION

Chromium oxide is reduced. $C_6H_{10} + C_7e^{6+} + H_2SO_4 \rightarrow C_7e^{3+}$

4. CALIBRATION OF THE TUBE

GAS CHROMATOGRAPHY

5. INTERFERENCE AND CROSS SENSITIVITY

Substance	Interference	Coexistence
Alcohols	Similar or brown stain is produced.	Higher readings are given.
Esters	"	"
Ketones	"	"
Aromatic hydrocarbons	"	"
Aliphatic hydrocarbons (more than C ₃)	Whole reagent is changed to pale brown.	If the maximum end point of the stain is discernable, the accuracy of readings is not affected.
Halogenated hydrocarbons	"	"